



Erratum: Vol. 51, No. RR-10

In the *MMWR Recommendations and Reports*, "[Guidelines for the Prevention of Intravascular Catheter-Related Infections](#)," published on August 9, 2002, on page 29 in Appendix B, an error occurred under the column heading "Replacement and relocation of device." The corrected table is as follows:

Appendix B

Summary of Recommended Frequency of Replacements for Catheters, Dressings, Administration Sets, and Fluids

Catheter	Replacement and relocation of device	Replacement of catheter site dressing	Replacement of administration sets	Hang time for parenteral fluids
Peripheral venous catheters	In adults, replace catheter and rotate site no more frequently than every 72–96 hours. Replace catheters inserted under emergency basis and insert a new catheter at a different site within 48 hours. In pediatric patients, do not replace peripheral catheters unless clinically indicated.	Replace dressing when the catheter is removed or replaced, or when the dressing becomes damp, loosened, or soiled. Replace dressings more frequently in diaphoretic patients. In patients who have large bulky dressings that prevent palpation or direct visualization of the catheter insertion site, remove the dressing and visually inspect the catheter at least daily and apply a new dressing.	Replace intravenous tubing, including add-on devices, no more frequently than at 72-hour intervals unless clinically indicated. Replace tubing used to administer blood, blood products, or lipid emulsions within 24 hours of initiating the infusion. No recommendation for replacement of tubing used for intermittent infusions. Consider short extension tubing connected to the catheter to be a portion of the device. Replace such extension tubing when the catheter is changed.	No recommendation for the hang time of intravenous fluids, including nonlipid-containing parenteral nutrition fluids. Complete infusion of lipid-containing parenteral nutrition fluids (e.g., 3-in-1 solutions) within 24 hours of hanging the fluid. Complete infusion of lipid emulsions alone within 12 hours of hanging the fluid. Complete infusions of blood products within 4 hours of hanging the product.
Midline catheters	No recommendation for the frequency of the catheter replacement.	As above.	As above.	As above.
Peripheral arterial catheters	In adults, do not replace catheters routinely to prevent catheter-related infection. In pediatric patients, no recommendation for the frequency of catheter replacement. Replace disposable or reusable transducers at 72-hour intervals. Replace continuous flush device at the time the transducer is replaced.	Replace dressing when the catheter is replaced, or when the dressing becomes damp, loosened, or soiled, or when inspection of the site is necessary.	Replace the intravenous tubing at the time the transducer is replaced (i.e., 72-hour intervals).	Replace the flush solution at the time the transducer is replaced (i.e., 72-hour intervals).
Central venous catheters including peripherally inserted central catheters and hemodialysis catheters	Do not routinely replace catheters.	Replace gauze dressings every 2 days and transparent dressings every 7 days on short-term catheters. Replace the dressing when the catheter is replaced, or when the dressing becomes damp, loosened, or soiled, or when inspection of the site is necessary.	Replace intravenous tubing and add-on devices no more frequently than at 72-hour intervals. Replace tubing used to administer blood products or lipid emulsions within 24 hours of initiating the infusion.	No recommendation for the hang time of intravenous fluids, including nonlipid-containing parenteral nutrition fluids. Complete infusions of lipid-containing fluids within 24 hours of hanging the fluid.
Pulmonary artery	Do not replace catheter to	As above	As above	As above

necessary.

Pulmonary artery catheters

Do not replace catheter to prevent catheter-related infection.

As above.

As above.

As above.

Umbilical catheters

Do not routinely replace catheters.

Not applicable.

Replace intravenous tubing and add-on devices no more frequently than at 72-hour intervals. Replace tubing used to administer blood products or lipid emulsions within 24 hours of initiating the infusion.

No recommendation for the hang time of intravenous fluids, including nonlipid-containing parenteral nutrition fluids. Complete infusion of lipid-containing fluids within 24 hours of hanging the fluid. Includes nontunneled catheters, tunneled catheters, and totally implanted devices.

[Return to top.](#)

Use of trade names and commercial sources is for identification only and does not imply endorsement by the U.S. Department of Health and Human Services.

References to non-CDC sites on the Internet are provided as a service to *MMWR* readers and do not constitute or imply endorsement of these organizations or their programs by CDC or the U.S. Department of Health and Human Services. CDC is not responsible for the content of pages found at these sites. URL addresses listed in *MMWR* were current as of the date of publication.

Disclaimer All *MMWR* HTML versions of articles are electronic conversions from ASCII text into HTML. This conversion may have resulted in character translation or format errors in the HTML version. Users should not rely on this HTML document, but are referred to the electronic PDF version and/or the original *MMWR* paper copy for the official text, figures, and tables. An original paper copy of this issue can be obtained from the Superintendent of Documents, U.S. Government Printing Office (GPO), Washington, DC 20402-9371; telephone: (202) 512-1800. Contact GPO for current prices.

**Questions or messages regarding errors in formatting should be addressed to mmwrq@cdc.gov.

Page converted: 8/15/2002

[HOME](#) | [ABOUT MMWR](#) | [MMWR SEARCH](#) | [DOWNLOADS](#) | [CONTACT](#)
[POLICY](#) | [DISCLAIMER](#) | [ACCESSIBILITY](#)

SAFER • HEALTHIER • PEOPLE™

Morbidity and Mortality Weekly Report
 Centers for Disease Control and Prevention
 1600 Clifton Rd, MailStop K-95, Atlanta, GA 30333,
 U.S.A



Department of Health
 and Human Services



This page last reviewed 8/15/2002